


Minta[®] Disposable RF Electrodes

Improves the Safety of the Patient and
The Medical Staff by Eliminating Cross
Contamination



Minta[®] RF Electrodes are used worldwide with Neurotherm, Radionics Baylis, Owl (Diros), Stryker, Cosman and other RF Lesion Generators previously supplied by and labelled as Neurotherm.

PAJUNK[®]

Minta Medical have joined together with **Pajunk[®]** the leading Cannula Manufacturer to offer Pain Management Clinics the Highest Quality RF Cannula and Electrode available.

Minta Medical

Leading the way in Pain Management Electrodes

CE 0086



MD81405

Minta® Disposable RF Electrodes and Pajunk® RF Cannula International Ordering Information

RF Lesion Generator	Minta® Disposable RF Electrode	Minta® Intermediate Cable	Pajunk® RF Cannula
Neurotherm JK2, 3, 4 & 25 NT1000 & NT1100	GABI - 50L GABI - 100L GABI - 150L	GABI-NS GABI-NS GABI-NS	001165 - 10B 001165 - 20B 001165 - 30BM
Radionics RFG - 3C or later	GABI - 50L GABI - 100L GABI - 150L	GABI-RS GABI-RS GABI-RS	001165 - 10B 001165 - 20B 001165 - 30BM
Baylis PMG - 115	GABI - 50L GABI - 100L GABI - 150L	GABI-BS GABI-BS GABI-BS	001165 - 10B 001165 - 20B 001165 - 30BM
Owl (Diros) URF - 3AP	SONNY - 50 SONNY - 100 SONNY - 150	Owl (Diros) C463 - 103 - BPTC -S	001165 - 10B 001165 - 20B 001165 - 30BM
Stryker	CLARI - 50 CLARI - 100 CLARI - 150	CLARI - ST CLARI - ST CLARI - ST	001165 - 10B 001165 - 20B 001165 - 30BM
Cosman RFG1A & 1B G4	GABI - 50L GABI - 100L GABI - 150L	GABI - CS GABI - CS GABI - CS	001165 - 10B 001165 - 20B 001165 - 30BM

PAJUNK also offer Special Active Tip Lengths and Curved / Angled RF Cannula

www.mintamedical.co.uk
www.pajunk.com