

# Minta® Disposable CORDOTOMY RF Electrode Kit

- Disposable electrodes improves the safety of the patient and medical staff by eliminating cross contamination, when performing CORDOTOMY procedures
- Precision tip for safe and easy insertion
- High accuracy calibration for consistent lesions
- High integrity electrode coating for maximum RF protection
- Radio opaque hub on needle for clearer x-ray visibility



**Minta® Disposable CORDOTOMY RF Electrodes** are used worldwide with Neurotherm, Radionics and other RF Lesion Generators.

**PAJUNK®**

Minta Medical have joined together with Pajunk® the leading Cannula Manufacturer to offer Pain Management Clinics the Highest Quality Lumbar Needle and Disposable CORDOTOMY Electrode available.

CE 0086



## Minta Medical

Leading the way in Pain Management Electrodes

## International Ordering Information for Minta® Disposable CORDOTOMY RF Electrode Kit

RF Lesion Generator	Minta® Disposable CORDOTOMY RF Electrode Kit	Minta® Intermediate Cable	Pajunk® Lumbar needle (included in kit)
Neurotherm JK2, 3, 4 & 25	MAX-T	MAX-N	1149-3Go80
Neurotherm NT1000 NT1100 NT2000	MAX-T	MAX-NT	1149-3Go80
Radionics RFG-3c plus RFG-3D	MAX-T	GABI-RS	1149-3Go80

Contact **Minta®** for details of using **Minta® CORDOTOMY Electrodes** with other RF machines

**“MANY THOUSANDS OF PATIENTS HAVE BEEN TREATED WORLDWIDE USING MINTA® CORDOTOMY RF ELECTRODES”**

**[www.mintamedical.co.uk](http://www.mintamedical.co.uk)**  
**[www.pajunk.com](http://www.pajunk.com)**